

## KTP 180X

The KTP 180X is a small, rotary pharmaceutical tablet press, delivering state-of-the-art performance for R&D applications, clinical batches and small scale production.

The machine is available in 2 distinctive configurations:

- KTP 180X Production
- KTP 180X Laboratory

### Highlights

- Interchangeable turret
- Mixed turrets for B+D tooling available
- Easy access to all product contact parts (weight of individual parts < 10kg)
- Unique scale-up functionality via scientifically designed approach
- Patented wear-free brake magnets

- 2 paddle force feed fill shoe available
- Individual tablet rejection system by compressed air
- Fully automated machine including a weight regulation
- Operation via a 12" integrated touch panel HMI
- R&D machine version monitors: upper and lower compression forces, ejection force, pre-compression force, tablet scraping force, upper and lower punch displacement
- Magnesium stearate lubrication available
- Data acquisition system integrated in the operation software
- Smallest footprint amongst R&D rotary tablet presses



**Technical data**

KTP 180X

Machine configurations	Production; Laboratory				
Number of press stations	8 (4+4)	8	13	16	19
Tool type (EU and TSM)	B+D	D	B+D	B	B
Die type	B+D	D	B+D	B	BB
Maximum tablet diameter (mm)	18 / 25	25	18 / 25	16 (18 for R&D)	13
Maximum die filling (mm)	17				
Maximum output* (tablets/hour)	16,800	33,000	54,600	67,200	79,800
Maximum pre-compression force (kN)	10				
Maximum main compression force (kN)	50**				
Main Drive Power (kW)	1.7				
Standard voltage (V), frequency (Hz)	400, 50/60				
Compressed air (bar)	4-6				
Weight (kg)	560				

\* The output depends on tablet size and compressed material

\*\* Maximum punch load dependant on tablet design

## Dimensions



