



### *Series LF 200 Romaco Macofar*

The new and innovative Romaco Macofar LF 2000 series is designed for high quality filling of pharmaceutical, nutraceutical and cosmetics liquid products into a wide range of bottles and vials, from 2 to 500 ml, at speeds of up to 12,000 units per hour.

#### **Highlights**

- Maintenance free
- Linear configuration
- Positive bottle transport system
- Bottle/vial bottom up filling
- Fast product change-over
- Customized line configurations
- High operational availability
- Excellent Overall Equipment Efficiency
- Low cost of ownership
- Very competitive price/ performance rating

<b>Technical Data</b>	<b>LF 200 FS</b>	<b>LF 200 FD</b>
Description	Filling + Single capping	Filling + Double capping
Motion	Continuous	Continuous
Max output (units/h)	12,000	12,000
Filling range (min-max), in ml	1 – 250	1 – 250
Bottle diameter (min-max) in mm	20 – 62	20 – 62
Bottle height, in mm	40 – 200	40 – 200
Installed power, in kW	12	18
Compressed air consumption, in NI/min	200	200
Compressed air pressure, in bar	6	6
Width, in mm	1,490	1,490
Length, in mm	4,990 (including capping module and conveyors)	5,510 (including capping module and conveyors )
Height, in mm	2,100	2,100
Weight, in kg	3,550	3,800

