

*The new generation of
automatic capping machines*

The capping machines in the Romaco Macofar MAC series integrate the most advanced technology for sealing cylindrical glass vials with aluminum flip-off, alu-caps. The machines comply with all established guidelines and regulations for aseptic processes.



Technical Data	MAC 1	MAC 2	MAC 4
Vials Diam., mm (min./max)	16 - 68		
Vials H., mm (min./max)	35 - 140		
Max. Output, vph	6.000	12.000	24.000
Closing heads	1	2	12
Motion	Intermittent	Intermittent	Continuous
Rated power, kW	2.5		
Weight, kg	800		

