



## 960 from Romaco Noack

The Noack solutions impress with excellent OEE values. Romaco continuously cuts customer's unit packaging costs by delivering optimal value for money. A balcony structure ensuring cGMP compliance together with integrated Track & Trace and containment solutions are the hallmarks of this technology. High level of automation with servo and auto controls ensures precise execution independent of the operator.

### Highlights

- QuickFeed – interchangeable feeding unit for product changeovers to be completed within seconds
- QuickAdjust – the automatic station control system for extremely short run-in times and highly reliable processes
- QuickChange – complete changeover in less than 30 minutes thanks to fail-safe control, easy plug and play, and ergonomic lightweight construction
- QuickTransfer – servo-driven blister transfer on a vacuum assisted overhead conveyor
- QuickClean – significantly shorter cleaning times due to easy access, reduced surfaces and a design without edges and corners
- QuickAdmin – absolutely reproducible processes and ergonomic HMI panel with intuitive visualisation and navigation
- QuickControl – maximum production safety along with Track & Trace, audit trail and batch records
- QuickConnect – Ready for Remote Assist; OPC UA interface for production data acquisition (PDA)

<b>Technical Data</b>	Noack 960
Output blister/minute max.	600
Blister tracks max.	6
Forming cycles/minute max.	100
Foil width max., in mm	235
Forming area max., in mm	225 x 180
Cutting area max., in mm	225 x 140
Forming depth max., in mm	12
Forming foil reel diameter max., in mm	700
Lidding foil reel diameter max., in mm	400
Machine dimensions (L x W x H), in mm	6,600 x 1,600 x 2,400
Weight, in kg	5,500
Remote Assist connection	Ethernet
PDA interface	OPC UA

Technical changes reserved

