



## FlexTrans from Romaco Siebler

The Siebler FlexTrans is the perfect solution for connecting primary packaging machines with a variety of different cartoners. With a 2 to 8 lane operation (depending on type and size of product) up to 900 packs per minute are handled. A wide range of pack sizes can be processed. The transfer/stacking technology allows fully automated parameter setting from the stored format recipe.

### Highlights

- Accurate Pick & Place operation
- Camera guided pick-up
- Rejection of non-qualified packs
- Minimised stress on finished packs
- Various stacking patterns available
- Quickest format changeover
- Individually pluggable suction heads
- Adjustable stacking unit

### Technical Data

	1-200	1-400	1-600	2-450	2-900
Single pack size min., in mm	40 x 50 x 3				
Single pack size max., in mm	120 x 240 x 13				
Speed, packs/minute max.	200	400	600	450	900
Remote Assist connection	Ethernet				
PDA interface	OPC UA				

