

### *HM 1 series from Romaco Siebler*

This latest generation strip packaging machine combines high performance "sleeved-shell" heating technology with its unique "QuickSeal" sealing technology and features quick release format and sealing modules to facilitate easy access for cleaning and quick product changeovers.

#### **Highlights**

- High availability
- Quick production start
- Modular format toolings
- Top sealing quality (repeatable/pre-set)
- Reproducible sealing parameters
- Perfect accessibility
- Small footprint
- Ready for Remote Assist
- OPC UA interface for production data acquisition (PDA)



<b>Technical Data</b>	HM 1-230 eco	HM 1-230	HM 1-350	HM 1-600
Foil width max., in mm	240	240	360	600
Product lanes max.	10	10	14	20
Output tablets/minute max.	3,000	3,500	4,900	7,000
Remote Assist connection	Ethernet			
PDA interface	OPC UA			

Technical changes reserved

