



E150 PLUS

The E150 PLUS is an economical and reliable, tablet press for various applications. It guarantees low service costs combined with an excellent tablet quality.

The machine is available in two distinctive configurations:

- E150 PLUS Pharma
- E150 PLUS Industrial - Food

Highlights

- Fully accessible change-point for quick tool change-over
- Availability of different fill shoe systems for powder feeding optimisation
- Automatic overload safety device – protects machine and tooling in case of over-pressure
- Patented wear-free punch brake magnets for constant redosing
- Integrated control cabinet and HMI
- Operation via 12" touch panel
- Monitoring of all process parameters

Technical data	E150 PLUS				
Machine configurations	Pharma; Industrial - Food				
Number of press stations	36	32	26	21	21
Tool type* (EU and TSM)	B	B	B	D	KILIAN D
Die type	BBS	BB	B	D	bounded die 38/23
Maximum tablet diameter (mm)	11	13	16	25	25
Maximum output (tablets/hour)**	225,500	200,400	162,800	131,500	131,500
Maximum die filling (mm)	15			20	45
Maximum pre-compression force (kN)	23				
Maximum main compression force (kN)	80***			100***	
Main Drive Power (kW)	4				
Standard voltage (V), frequency (Hz)	400, 50/60				
Compressed air (bar)	6				
Weight (kg)	1,600				

* Special tools/formats available on request

** The output depends on tablet size and compressed material

*** Maximum punch load dependant on tablet design

Dimensions



