

HM 1 series from Romaco Siebler

This latest generation strip packaging machine combines high performance “sleeved-shell” heating technology with its unique “QuickSeal” sealing technology and features quick release format and sealing modules to facilitate easy access for cleaning and quick product changeovers.

Highlights

- High availability
- Quick production start
- Modular format toolings
- Top sealing quality (reproducible setting parameters)
- Ready for Industrie 4.0
- Perfect accessibility
- Small footprint
- Ready for Remote Assist
- OPC UA interface for production data acquisition (PDA)



Technical Data	HM 1-200	HM 1-230 eco	HM 1-230	HM 1-350	HM 1-600
Maximum foil width (mm)	200	240	240	360	600
Maximum no.of lanes	4	10	10	15	20
Maximum output (tablets/minute)	1,200	3,000	3,500	5,250	7,000
Maximum output (capsules/minute)	800	1,600	1,800	3,000	4,000
Remote Assist connection	Ethernet				
PDA interface	OPC UA				

Technical changes reserved

