

Capsule Filling Machine CD 40 from Romaco Macofar

Romaco Macofar's compact, intermittentmotion capsule filling and closing machine is chosen by users all over the world for its flexibility and user-friendliness. The machine is ideal for medium production runs, with mechanical speed of up to 40,000 capsules per hour.

The dosators are vacuum assisted, which allows powders with difficult flow characteristics to be handled, giving improved compaction of powder.

Highlights

- Outstanding filling accuracy via verticalaction vacuum dosators
- Pellet dosage by means of a gravity system, preserving the pellet integrity and characteristics
- Capability to fill a wide range of powders and/or pellets
- Minimum format parts for fast changeover
- Statistical weight control system for filled capsules
- Optional automatic feeding groups for products and capsules
- Full technical documentation for validation purposes



Capsule size	00 - 0 - 1 - 2 - 3 - 4 - 5	
Rated power (kW)	3.5	
Machine compressed air (NI/min)	150	
Machine + CW compressed air (NI/min)	350	
Machine dimensions, approx. (W \times H \times L, in mm)	980 x 1.510 x 2.080	
Weight, approx. (kg)	1,200	

CD 40

40,000

Technical Data

Maximum output (capsules/hour)

