

*The new generation of
automatic capping machines*

The capping machines in the Romaco Macofar MAC series integrate the most advanced technology for sealing cylindrical glass vials with aluminum flip-off, alu-caps. The machines comply with all established guidelines and regulations for aseptic processes.



Technical Data	MAC 2	MAC 4
Vials diameter (mm), min.–max.	16–35	
Vials height (mm), min.–max.	35–100	
Speed up to (vials/hour)*	12,000	24,000
Closing heads	2	12
Motion	Intermittent	Continuous
Rated power (kW)	2.5	
Weight (kg)	800	

* Depending on product and machine configuration

